








PROSTHETIC COMPONENTS - KNEE JOINTS

	ORDER NO.	DESCRIPTION
	KN-1	Single Bar Pneumatic Knee Unit w/ Stance Control
	KN-3	4 Bar mechanical Knee Unit w/ carbon graphite frame
	KN-4	4 Bar mechanical Knee Unit, aluminum frame, short
	KN-5	4 Bar Pneumatic Knee Unit w/ carbon graphite frame
	KN-8	4 Bar mechanical Knee Unit, aluminum frame, tall, KnDa
	KN-55	Titanium Weight Activated Brake Knee w/ extension assist(120 gm)
	MT99	AL 4 Hole Aluminum Lamination Ring / Screws

Toll Free - 1 (877) 395-0081

MOT Myrdal
Orthopedic
Technologies
Technical O&P Solutions



PROSTHETIC COMPONENTS - TITANIUM

	ORDER NO.	DESCRIPTION
	MT-10	Titanium 4 hole connector, pyramid 136 kg/300 lbs. Part Weight = 46 gm
	MT-11	Titanium 30mm tube clamp w/screws 136 kg/300 lbs. Part Weight = 79 gm
	MT-12	Titanium SACH foot adapter w/ bolt and lock washer 136 kg/300 lbs. Part Weight = 94 gm
	MT-13	Titanium 4 ear connector w/ pyramid 136 kg/300 lbs. Part Weight = 56 gm
	MT-14	Titanium Tube clamp assembly 400mm, receiver w/ screws 136 kg/300 lbs. Part Weight = 299 gm
	MT-15	Titanium Rotating pyramid receiver w/ screws 136 kg/300 lbs. Part Weight = 46 gm
	MT-16	Titanium Rotating pyramid adapter 136 kg/300 lbs. Part Weight = 50.5 gm






Toll Free - 1 (877) 395-0081

MOT Myrdal
Orthopedic
Technologies
Technical O&P Solutions

PROSTHETIC COMPONENTS - TITANIUM

ORDER NO.	DESCRIPTION
	<p>MT-17 Titanium 4 hole rotating adapter base w/ screws 136 kg/300 lbs. Part Weight = 70 gm</p>
	<p>MT-18 Titanium 4 hole connector w/ receiver and screws 136 kg/300 lbs. Part Weight = 67 gm</p>



PROSTHETIC COMPONENTS - STAINLESS STEEL

	<p>MT-71 SS. 30mm tube clamp w/ screws 114 kg/250 lbs. Part Weight = 128 gm</p>
	<p>MT-72 SS. SACH foot adapter w/ bolt and lock washer 114 kg/250 lbs. Part Weight = 110 gm</p>
	<p>MT-73 SS. 4 ear connector w/ pyramid 114 kg/250 lbs. Part Weight = 95 gm</p>
	<p>MT-74 SS. Tube clamp assembly 400mm, receiver w/ screws 114 kg/250 lbs. Part Weight = 327 gm</p>
	<p>MT-75 SS. Rotating pyramid receiver w/ screws 136 kg/300 lbs. Part Weight = 69.75 gm</p>






Toll Free - 1 (877) 395-0081

MOT Myrdal
 Orthopedic
 Technologies
Technical O&P Solutions



PROSTHETIC COMPONENTS - STAINLESS STEEL

ORDER NO.	DESCRIPTION
	<p>MT-76 SS. 3 prong rotating adapter base w/ screw 136 kg/300 lbs. Part Weight = 114 gm</p>
	<p>MT-80 SS. 4 hole connector w/ pyramid 114 kg/250 lbs. Part Weight = 104.5 gm</p>




PROSTHETIC COMPONENTS - ALUMINUM

	<p>MT-91 ALUMINUM 30mm tube clamp w/ screws 100 kg/220 lbs. Part Weight = 86.5 gm</p>
	<p>MT-92 ALUMINUM SACH foot adapter w/ bolt and lock washer 100 kg/220 lbs. Part Weight = 104 gm</p>
	<p>MT-93 ALUMINUM 4 ear connector w/ pyramid 100 kg/220 lbs. Part Weight = 64 gm</p>
	<p>MT-94 ALUMINUM Tube clamp assembly 400mm, receiver w/ screws 100 kg/220 lbs. Part Weight = 308 gm</p>
	<p>MT-100 ALUMINUM 4 hole connector w/ pyramid 100 kg/220 lbs. Part Weight = 63.5 gm</p>

PROSTHETIC COMPONENTS - FEET

	ORDER NO.	DESCRIPTION
	FT-1	SACH FOOT split toe Left or Right, sizes 24 - 28, 3/8" heel ht.
	FT-2	Ti. Single axis ankle unit w/ foot, sizes 24 - 28, 3/8" heel ht.
	FT-3	Dynamic energy storing foot w/ pyramid, sizes 23-28
	PY-1	Carbon graphite pylon, 30mm. 12 inches 136 kg/300 lbs. Part Weight = 86 gm

ORTHOTIC COMPONENTS

	LTSY01JT	NIAGARA Knee Joint - Left
	RTSY01JT	NIAGARA Knee Joint - Right Designed for single upright lateral placement using Hollister Alignment
	OEJ10520	Elbow Joint ROM Bent/Double Slotted Uprights, Black Anodized
	OEJ10520-C	Elbow Joint ROM Bent/Double Slotted Uprights, Clear Anodized
	OEJ105201	Elbow Joint ROM Bent/Double Slotted Uprights, Black/Hardware
	OEJ105201-C	Elbow Joint ROM Bent/Double Slotted Uprights, Clear / Hardware
	OEJ10522	Elbow Joint ROM - Double Slotted Uprights, Black Anodized
	OEJ10522-C	Elbow Joint ROM - Double Slotted Uprights, Clear Anodized
	OEJ10523	Elbow Joint ROM - Solid upright - no slots, Black Anodized
	OEJ10523-C	Elbow Joint ROM - Solid upright - no slots, Clear Anodized
	OEJFMPCK (image not shown)	Elbow Joint FM - Polycentric Solid upright Black Anodized Polycentric

Toll Free - 1 (877) 395-0081

MOT Myrdal
Orthopedic
Technologies
Technical O&P Solutions